








Biocard[®] HIV 1/2 RAPID TEST KIT

Biocard[®] HIV 1/2 Rapid test kit is an in-vitro immunochromatographic immunoassay for the qualitative simultaneous detection and differentiation of HIV-1 and HIV-2 antibodies(IgG, IgM and IgA)


-  Approved by National Institute Biologicals (NIB)
-  Sensitivity : 100% Specificity : 100%
-  Single test device to screen and differentiate HIV-1 & HIV-2
-  Sample types: Serum, Plasma (10µL)/ Whole Blood (20µL)
-  Precoated with gp41 & gp120 HIV-1 antigens and gp36 HIV-2 antigens
-  Results with in 20 minutes
-  Storage: 4°C to 30°C



Biocard®


HIV 1/2 RAPID TEST KIT

Test Procedure



Result Interpretation

For Biocard® HIV Device



Performance

Status	Positive	Negative
Positive	110	0
Negative	0	650
Total	110	650

Diagnostic Sensitivity: 100 %

Diagnostic Specificity: 100 %

Ordering Information

Product Code	Product Name	Pack Size	Storage	Shelf Life
30201136	Biocard® HIV 1/2 RAPID TEST KIT	30 Test	4-30°C	24 Months

*Please note: Product specifications are subject to change without prior notice owing to product modifications, improvements / up-gradation.



corporate@trivitron.com | www.trivitron.com

Customer Care Helpline: +91 98400 80008

