

EARLY DETECTION SPEEDS UP RECOVERY!



RNAsure RNA EXTRACTION KIT

- The Kit is based on Magnetic Bead Technology for extraction of RNA
- Manual Extraction Kit
- Storage at room temperature (15-25°C)

RNAsure RNA Extraction Kit

The RNAsure Corona virus (SARS-CoV-2) RNA extraction kit is designed for the isolation of RNA from nasal/throat swabs collected in viral transfer medium (VTM). This kit provides reagents and magnetic beads for isolation of 500 samples.


Specification

- RNAsure Kit is designed for rapid manual preparation of RNA from nasal/throat swabs collected in VTM
- The kit is designed for use with RNAsure Pure magnetic separator which can hold 12 or 16 micro-centrifuge tubes for the easy isolation of 12 or 16 RNA samples in less than 20 minutes
- Kit can be easily adoptable in any 96-well magnetic separator plate or other automated magnetic separation systems
- Storage - room temperature (15-25°C)
- Validity Period – Stable for 12 months from the date of manufacture

Kit Contents

Kit Content	For 500 Preparations
RNAsure Pure Beads	20ml
Lysis Buffer	50ml
Wash Buffer Concentrate	2x50ml
Elution Buffer	50ml
Carrier RNA	1 Vial
User Manual	1
Magnetic Separator (12 / 16 tubes)	1

Schematic Diagram of Process Flow



*Please note: Product specifications are subject to change without prior notice owing to product modifications, improvements / up-gradation.