










Biocard®


Dengue Ag (NS1) Rapid Test Kit

Biocard® Dengue Ag (NS1) Rapid Test Kit is a rapid immunochromatographic assay for the qualitative detection of Dengue NS1 Antigen in the human serum / plasma.


-  Technique: Rapid test based on lateral flow
-  Accuracy: Detection of Dengue in the seronegative stage i.e. from Day 1 of fever
-  Detection: All 4 Dengue serotypes (DEN-1, DEN-2, DEN-3 and DEN-4)
-  Results: 20 minutes
-  Longer shelf life: 24 months
-  Storage: 4°C to 30°C
-  Convenient pack size: 30 tests
-  Sensitivity: 99.06%
-  Specificity: 100%

Dengue Ag (NS1) Rapid Test Kit

Test Procedure



Results Interpretation



Performance

Status	Positive	Negative	Total
Positive	213	2	215
Negative	0	250	250
Total			465

Diagnostic Sensitivity: 99.06%

Diagnostic Specificity: 100%

Ordering Information

Product Code	Product Description	Pack Size	Specimen
30201448	Biocard® Dengue Ag (NS1) Rapid Test Kit	30 Tests	Serum/Plasma

*Please note: Product specifications are subject to change without prior notice owing to product modifications, improvements / up-gradation.



corporate@trivitron.com | www.trivitron.com

Customer Care Helpline: +91 98400 80008

