



# Personal Protective Equipment

Safety Solutions for Healthcare Professionals



**Elastic Hood**

Provides optimal protection & a comfortable fit

**Seam Sealed**

Prevents any sort of liquid penetration to maintain complete protection

**Longer Zipper**

Allows for easy donning and doffing

**Ample Crotch**

**Flap Cover**

**Elastic Wrist**

To enable proper range of movement

**Disposable**

**Kiran PPE Kit Includes:**

- Protective Coverall - 1 No | Shoe Cover - 1 Pair
- Hand Gloves - 1 Pair | Face Shield - 1 No
- Face Mask - 1 No | Disposable Bag - 1 No

**Resistance to Blood Penetration**

Complies to ASTM F 1670

**Standard**

**Optional**



**Kiran Face Shield Plus**

500 µm PET Clear Flexi sheet with Anti-Scratch Coating (Reusable)



**Face Mask**

3 Ply Face Mask Comfortable & Breathable (Disposable)



**Shoe Cover**

Same material as Coverall (Disposable)



**Gloves**

Nitrile Rubber, Latex-free, odour-free & powder-free (Disposable)



**Kiran Face Shield Pro**

2 mm thick Acrylic Visor with Anti Scratch Coating (Reusable)



**N95 Mask**

Breathable & Washable (Reusable)



**Liquid Resistant | Bound Seams | Soft Fabric**

### Kiran PPE Pro

#### Coverall Fabric Type

Microporous Polypropylene  
Non-woven Fabric with  
Extrusion Lamination (PP+PE)  
The coverall acts as a  
barrier to water, vapour, etc  
(Disposable)

### Kiran PPE Max

#### Coverall Fabric Type

Woven Polyester Fabric  
with PU Coating which  
creates a waterproof barrier  
(Disposable)

### Kiran PPE Plus

#### Coverall Fabric Type

Non-woven Fabric  
with TPU film lamination  
(Disposable)

\*Please note: Product specifications are subject to change without prior notice owing to product modifications, improvements / up-gradation.

# KIRAN

D -117, 134, T.T.C.Industrial Area  
Nerul, Navi Mumbai- 400706. INDIA

Ph: +91 22 27630942, Fax: +91 22 27671790

Email id : kiran@kiranxray.com, leads@trivitron.com

A Division of TRIVITRON Healthcare



[www.kiranxray.com](http://www.kiranxray.com)