



Personal Protective Equipment

Safety Solutions for Healthcare Professionals



Elastic Hood

Provides optimal protection & a comfortable fit

Seam Sealed

Prevents any sort of liquid penetration to maintain complete protection

Longer Zipper

Allows for easy donning and doffing

Ample Crotch

Flap Cover

Elastic Wrist

To enable proper range of movement

Disposable

Kiran PPE Kit Includes:

- Protective Coverall - 1 No | Shoe Cover - 1 Pair
- Hand Gloves - 1 Pair | Face Shield - 1 No
- Face Mask - 1 No | Disposable Bag - 1 No

Resistance to Blood Penetration

Complies to ASTM 1670

Standard

Optional



Kiran Face Shield Plus

500 µm PET Clear Flexi sheet with Anti-Scratch Coating (Reusable)



Face Mask

3 Ply Face Mask Comfortable & Breathable (Disposable)



Shoe Cover

Same material as Coverall (Disposable)



Gloves

Nitrile Rubber, Latex-free, odour-free & powder-free (Disposable)



Kiran Face Shield Pro

2 mm thick Acrylic Visor with Anti Scratch Coating (Reusable)



N95 Mask

Breathable & Washable (Reusable)



Liquid Resistant | Bound Seams | Soft Fabric

Kiran PPE Pro

Coverall Fabric Type

Microporous Polypropylene
Non-woven Fabric with
Extrusion Lamination (PP+PE)
The coverall acts as a
barrier to water, vapour, etc
(Disposable)

Kiran PPE Max

Coverall Fabric Type

Woven Nylon Fabric
with PU Coating which
creates a waterproof barrier
(Disposable)

Kiran PPE Plus

Coverall Fabric Type

Non-woven Fabric
with TPU film lamination
(Disposable)

*Please note: Product specifications are subject to change without prior notice owing to product modifications, improvements / up-gradation.

KIRAN

D -117, 134, T.T.C.Industrial Area
Nerul, Navi Mumbai- 400706. INDIA

Ph: +91 22 27630942, Fax: +91 22 27671790

Email id : kiran@kiranxray.com, leads@trivitron.com

A Division of TRIVITRON Healthcare



www.kiranxray.com