

SonoRad V8

Premium

Slim - Smart - Significant - Streamlined



Introducing

SonoRad V8

Premium

SonoRad V8 Premium, a slim console color doppler system offers ideal patient care with its easy mobility, smart features, significant performance and streamlined workflow to strengthen the clinical confidence and speed the diagnostic decision.


Excellent Design



- ▶ 19" LED Monitor
 - 90°~ 90° rotatable
 - 0°~ 90° foldable
- ▶ Replaceable probe holder
- ▶ Cable management
- ▶ Facing front printer socket
- ▶ Rotatable keyboard: -45°~ 45°
- ▶ Three probe connectors

Advanced Features

Demo Features




Complete DICOM function

DICOM SR, transfer, print, etc. Professional DICOM structure report; User definable measurement in DICOM SR; User Selectable DICOM transfer format;




HIP graf & image rotation



- ▶ A graph for hip orthotics diagnosis helps the doctor to make an easier and more accurate diagnosis during the pediatric hip scanning.
- ▶ Different angles indicate different levels of hip deformity, which becomes easier and more obvious to see with the aid of the graph (I, II, D, IIIa, IIIb).
- ▶ Image rotation: 0°, 90°, 180°, 270° for right direction of the HIP scanning and measurement



Auto IMT



- ▶ The automatic intima-media thickness measurement for fast carotid exam
- ▶ The detailed result of the intima thickness provides automatic diagnosis

Streamline Workflow

- ▶ Streamlined process, six in one key
- ▶ Shortcut key for quick reaction
- ▶ User -friendly interface
- ▶ Anatomy-specific presets
- ▶ Efficient data management



Wide Range of Probes



2.0 - 6.8 MHz
Convex



4.0 - 15.0 MHz
Linear



4.0 - 12.0 MHz
Transvaginal



4.0 - 15.0 MHz
Transvaginal




2.0 - 6.8 MHz
Micro-Convex




4.0 - 12.0 MHz
Micro-Convex


Clinical Images




Carotid, B Mode




Carotid, PW Mode




Finger Vessel, CPA Mode




Vertebral Artery, C Mode




Thyroid, B Mode




Liver, C Mode



Fetal Heart, B Mode



Umbilical Cord, C Mode



Kidney, C Mode

TRIVITRON

HEALTHCARE

speaking your language

www.trivitron.com
corporate@trivitron.com | +91 98400 80008

