

Storage Cabinet

Technical Specification for Tissue-Tek® Ultra™ II Cabinet (Item Code: 3990)

- Dimensions (WXDXH) 59.5 X 71.3 X 110 cm
- Weight
 - Empty 123 Kg
 - Filled with Slides 585 Kg
 - Filled with Blocks 310 Kg
- Sturdy steel construction body ensures long-term use
- Drawer lock prevents possible tampering
- Each drawer has a plastic covered labeling tab for positive identification
- Convenient removable trays accommodate upto 96000 slides or 22000 blocks
- Coding strips can be used to separate slides or blocks by numbers
- Filing cabinets can be bolted together and stacked two high for maximum storage capacity