



The Compass™ Stainer

Optimized Staining Throughput with
an Efficient Laboratory Footprint

Compass™
STAINER



HOLOGIC® |  **GEN-PROBE**

The Compass™ Stainer

Before the introduction of the Compass™ Stainer, ThinPrep® with Imaging for small and mid-sized laboratories required a large footprint staining device that made inefficient use of valuable bench space. Laboratory managers were forced to make difficult trade-offs when allocating resources.

Hologic® solved this issue by developing the Compass Stainer to **enable laboratories of all sizes** to choose ThinPrep with Imaging for a technological advancement in the detection of cervical disease.

The Compass Stainer is **CE marked**, and is **validated** for the ThinPrep® Imager stain protocol.

A fully automated stainer for histologic and cytologic slides, the Compass Stainer provides users with:

- A **compact footprint** to maximize bench space
- The ability to run 20 or 30 position slide racks for **flexible, LEAN processing**
- Up to four stain racks can run simultaneously
- **Efficient use of consumables**, each reagent bath holds 350 to 360 ml so stain isn't wasted
- A **flexible programming model** supports a wide range of applications



Processed Slides	Time Period
30 Slides	28 minutes
120 Slides	85 minutes
750 Slides	1 work day

The complete solution for small and mid-sized laboratories



Compass Stainer Specifications¹

Instrument Dimensions	Width: 81 cm/32", Height: 41 cm/16", Depth: 79 cm/31"
Weight	25 kg / 55 lbs
Operating Temperature	15 – 35°C / 59 – 95°F
Operating Humidity	20% – 80% RH non-condensing
Electrical Voltage	220/240 VAC, 1 amp rated outlet 100/120 VAC
Power/Frequency	Maximum 180 watts at 50 – 60 Hz

¹Compass Stainer Operator's Manual, Part Number MAN-03101-001.

At Hologic, we turn passion into action, and action into change.

Hologic is defining the standard of care in women's health. Our technologies help doctors see better, know sooner, reach further and touch more lives.

BREAST IMAGING SOLUTIONS • INTERVENTIONAL BREAST SOLUTIONS • BONE HEALTH
PRENATAL HEALTH • GYNECOLOGIC HEALTH • MOLECULAR DIAGNOSTICS



www.hologic.com

North America / Latin America

35 Crosby Drive
Bedford, MA 01730-1401
USA

Europe

Otto-von-Guericke-Ring 15
65205 Wiesbaden-Nordenstadt
Germany

Asia Pacific

7th Floor, Biotech Centre 2
No. 11 Science Park West
Avenue
Hong Kong Science Park
Shatin, New Territories
Hong Kong

Australia / New Zealand

Suite 402, Level 4
2 Lyon Park Road
Macquarie Park NSW 2113
Australia