

iBeat3V Pro



ECG

- Simultaneous 12 leads ECG acquisition and display
- Simultaneous 3 leads ECG printing with 80mm thermal paper
- 7 inch 800 x 480 TFT color LCD
- Alpha-numeric keyboard
- Built-in rechargeable lithium battery
- Built-in 4G SD card, save up to 10,000 ECG reports in XML/PDF/BMP format (optional)
- Pacemaker detection
- Precise ECG auto-interpretation program
- USB SLAVE port can be connected to PC-ECG software (optional)
- USB MASTER port can be connected to an external USB printer supporting PCL6 language to print 12 leads ECG reports on A4 paper (optional)
- Network transmission (optional)
- Software upgrade automatically via SD card

SPECIFICATIONS*

The recording paper with 80mm width, the clearly printed reports include ECG measurement data, patient ID, sex, age etc

3CHx4



3CHx4+1R



- Lead** : Standard 12 leads
- Input Circuit** : Built-in defibrillator protection
- Acquisition Mode** : Simultaneous 12 leads acquisition
- Measurement Range** : $\pm 5\text{mVpp}$
- Time Constant** : $\geq 3.2\text{s}$
- Frequency Response** : 0.05Hz~150Hz (-3dB)
- Sensitivity** : 2.5 mm/mV, 5 mm/mV, 10 mm/mV, 20mm/mV
- Input Impedance** : $\geq 50\text{M}\Omega$
- Input Circuit Current** : $< 0.05\mu\text{A}$
- Noise Level** : $< 15\mu\text{Vp-p}$
- Anti Baseline Drift** : Automatic
- Patient Current Leakage** : $< 10\mu\text{A}$ (220V-240V)
- Filter** : EMG Filter: 30 Hz/40Hz
DFT Filter : 0.05 Hz/0.15 Hz/0.25 Hz/0.5Hz
AC Filter : 50Hz/60Hz



Main Interface



File Management Interface

- CMRR** : 98dB
- Recording Mode** : Automatic: 3CHx4+1R, 3CHx4
Manual: 3CH, 1CH, 1CH+1R
Rhythm: Any Lead selectable
- Recording System** : Thermal matrix printing
- Recording Paper** : Roll: 80mmx20m
- Paper Speed** : 6.25 mm/s, 12.5 mm/s, 50mm/s
- Display** : 800x480 pixel TFT LCD
- Interface** : Network interface, USB Master/Slave, SD card
- Safety Standard** : IEC I/CF
- Sample Rate** : Normal : 1000sps/channel
Pacemaker : 10,000sps/channel
- Pacemaker Detection** : $\geq \pm 2\text{mV}$ $\geq \pm 0.1\text{ms}$
- Dimension** : 310mmx244mmx65mm
- Weight** : Approx. 2.25Kg (Include battery)
- Power Supply** : AC : 100~240V, 50/60Hz, 70VA~100VA
DC : Built in Rechargeable Li-ion battery



- Head Office** : # 15, IV Street, Abhiramapuram, Chennai -18
Ph : +91 44 24985050 (6 lines) Fax : +91 44 24985757 / 24672782
- Regional Offices** : BANGALORE (+91 8022115491) | DELHI (+91 11 41354650)
KOLKATA (+91 33 24224493) | MUMBAI (+91 22 66612526 / 27 / 28)
- Branch Offices** : Chandigarh, Cochin, Coimbatore, Jaipur, Nagpur, Pondicherry & Pune
- International Offices** : UAE, Oman, Sri Lanka, Seychelles, Maldives, Mauritius, Bangladesh & Nepal

TRIVITRON
HEALTH CARE
speaking your language

Certified ISO 9001:2004 Company

www.trivitron.com

+91 - 98400 80008 (24 Hour Helpline)

*Please note: Product specifications are subject to change without prior notice owing to product modifications, improvements / up-gradation.