

URIT *Since 1984*

Global Diagnostics Supplier



Product Catalogue

URIT Medical Electronic Co., Ltd.

Head Office

Add. Add.: No. D-07 Information Industry District, High-tech Zone, Guilin, Guangxi 541004, P. R. China
Tel: +86-773-2288586 Fax: +86-773-2288560

Shenzhen Office

Add.: Room EJ, 24th Floor, New Baohui Building, Nanhai Road, Nanshan District, Shenzhen, P.R. China
Tel: +86-755-26050853 Fax: +86-755-86036150
E-mail: export@uritest.com service@uritest.com [Http://www.urit.com](http://www.urit.com)

Version:20170710



01/02
Company
Profile

03/06
Hematology
Analyzers
URIT-5500
URIT-5250
URIT-5380
URIT-5160
URIT-3000Plus
URIT-2900Plus

07/10
Urine
Analyzers
URIT-1600
URIT-1280
URIT-1000Plus
URIT-560
URIT-500B/500C
URIT-180
URIT-50

11/13
Chemistry
Analyzers
URIT-8280
URIT-8210
URIT-880
URIT-810

14/16
Point-of-care
Testing
URIT-12
URIT-15
URIT-10
URIT-28
URIT-25
URIT-31

17/17
Electrolyte
Products
URIT-910C Plus

17/17
Semi-c
coagulation
Products
URIT-600/610

18/18
Microplate
Reader
& Washer
Products
URIT-660
URIT-670

19/20
Medical
Software
LIS
Laboratory
Information
System

21/22
Reagent
Hematology Reagent
& Control
Chemistry Reagent
& Control
Urine Reagent Strips
UQ Series Control
Electrolyte Reagent

23/26
Veterinary
Product
URIT-5250 Vet
URIT-3000Vet Plus
URIT-2900Vet Plus
URIT-50 Vet
URIT-31Vet
URIT 11V Vet
URIT-880 Vet
URIT-9000 Vet

27/28
Medical
Electronic
Product
URIT-9000
URIT-7200



Company Profile

URIT is a global diagnostics supplier, with the mission to significantly improve the quality and efficiency for laboratories and hospitals around the world.

Since established in 1984, URIT has medical equipment production base of 80000 square meters and a number of imported production lines for manufacturing advanced automated machine and reagent; it is a manufacturer of China's high-end medical equipment research and development.

URIT's product range covers hematology analyzers, chemistry analyzers, POCT, urine analyzers, reagents, medical software, veterinary products and medical electronic products.

R&D Team & Achievements

Intelligent R&D team with over 200 outstanding engineers and specialists

Obtaining core technologies & patents on all of URIT's products

Quality Commitment

Strictly conducted in compliance with ISO9001 and ISO13485 high-quality products with CE & FDA approval

Marketing Network

Domestic market: More than 70 branch offices, holding leading market share

Overseas market: Marketing prosperously in more than 160 countries


Service Guarantee

Over 400 appointed domestic and 70 appointed overseas service stations

Professional service team, immediate response service and effective products training

Inspired by the vision of To be the World's Leading Diagnostics Supplier, Providing Excellent Services!, URIT is becoming a leading enterprise for diagnostics.

URIT devotes to presenting Chinese leading technology to the world. We warmly welcome friends from around the world cooperating together with URIT to bring international advanced technology back to China.



Hematology Analyzers

URIT-5500

5-Diff Hematology Analyzer

- 28 parameters including RET%, RET#, IRF
- Up to 110T/H
- Automatic sampling with a capacity of 120 samples
- Laser light multi-dimensional cell classification & flow cytometry technology
- Double mode for WBC counting (Impedance Count and Optical Count)
- High stability, long using life He-Ne Lasers
- Ceramic sample rotatory valve(SRV)quantified system with high accuracy and efficiency
- 3D Image Analysis
- WBC 5-diff with sheath reagent only
- Whole computer data station and display
- Large storage capacity:200,000 samples
- Supprt LIS and HIS with HLT protocol
- Internal barcode reader



URIT-5250

5-Diff Hematology Analyzer

- 28 parameters including RET%,RET#,IRF
- Up to 60T/H
- Laser light multi-dimensional cell classification & flow cytometry technology
- Doule mode for WBC counting
- Resistant RBC mode
- WBC 5-diff with sheath reagent only
- Support both whole blood and capillary blood samples with 20ul volume



URIT-5380

5-Diff Hematology Analyzer

- 28 parameters including RET%, RET#, IRF
- Laser light multi-dimensional cell classification
- Double mode for WBC counting (Impedance Count and Optical Count)
- Resistant RBC mode
- Up to 60T/H
- Internal barcode reader



URIT-3000Plus

Automatic Hematology Analyzer

- 60T/H, 21 parameters and 3 histograms
- 10.4 inch color LCD display & Linux operation system
- More than 200,000 data storage capacity
- Excellent data management
- Support USB data backup & system upgrade
- High efficient self checking system, less maintenance
- Support LIS and HIS with HI7 protocol
- Touch screen optional



URIT-5160

5-Diff Hematology Analyzer

- 24 parameters
- Laser light multi-dimensional cell classification
- Double mode for WBC counting (Impedance Count and Optical Count)
- Resistant RBC mode
- Up to 60T/H
- WBC 5 part diff with sheath reagent only
- Touch screen



URIT-2900Plus

Automatic Hematology Analyzer

- 30T/H, 21 parameters and 3 histograms
- 10.4 inch color LCD display & Linux operation system
- More than 200,000 data storage capacity
- Excellent data management
- Support USB data backup & system upgrade
- Support LIS and HIS with HL7 protocol
- Touch screen optional





URIT Fully Automatic Urinalysis System

- With 32 years experience in urinalysis product development
- Fully automatic urinalysis system consists of URIT-1600, URIT-1280, connecting bridge, work station
- Proficiency and efficiency in urinalysis in the most convenient way
- Adopting the most advanced image capture and multi-core parallel processing technique

Urine Analyzers

URIT-1600 Automatic Urine Analyzer

- Up to 300T/H
- Up to 14 parameters including Cr, MA, Ca, VitC
- 8 inch color display with touch screen
- Integrated barcode reader for sample
- More than 100,000 data memory capacity
- Accurate sample quantitative drop system
- Stable strip reaction time control
- COM and LAN port for LIS and HIS communication



URIT-1280 Automatic Urine Sediment Analyzer

- Three counting channels in flow cell
- Whole view filed microscopic images of sediment with copmuter
- 11 parameters: RBC,WBC,SQEP,NSE,HYA, UNCC,CAOX,BACT,YST, MUCUS,URIC, etc.
- Up to 100T/H and HIS communication
- Walk-away automation, high-efficiency washing system



URIT-1000Plus Automatic Urine Sediment Analyzer

- Available for preset the 9 test items
- Phase detection of RBC in urine
- Flat-field achromatic automatic continuous zoom microscopic test system
- Direct sampling with urine cup
- Sample volume only requires 1.2 ml
- Double channels to test



URIT-500B/500C Urine Analyzer

- Touch screen(URIT-500B) / Keypad(URIT-500C), FDA approved (URIT-500B)
- Excellent throughput, up to 500T/H
- Up to 14 parameters including Cr, MA, Ca, VitC
- High automation, continuous urine strip loading
- Easy maintenance design for daily cleaning
- Barcode reader optional
- URIT UQ-10, UQ-11 and UQ-14 control available



URIT-560 Urine Analyzer

- Up to 600T/H
- Up to 14 parameters including Cr, MA, Ca, VitC
- CIS detectors, faster and more accurate results
- Input more patient information
- Barcode reader optional, touch screen
- URIT UQ-10, UQ-11 and UQ-14 control available



URIT-180 Urine Analyzer

- Up to 120T/H and 11 parameters available
- User-friendly operation
- Friendly operation with large LCD screen
- Barcode reader optional
- URIT UQ-10 & UQ-11 UQ-14 control available



URIT-50 Urine Analyzer

- Up to 120T/H and 11 parameters available
- FDA approved
- Portable
- Barcode reader optional
- URIT UQ-10 & UQ-11 UQ-14 control available





Chemistry Analyzers

URIT-8280

Automatic Chemistry Analyzer

- Random Access, constant 640T/H, 1000T/H with ISE
- High accuracy grating optical system
- Durable ceramic syringes, 4 unattached probes, 2 high-performance mixers
- 24 hours non-stop water cycle cooling system to keep the reagent at 2-8°C
- Metal solid heating for the 37°C incubation
- Support LIS interface
- Barcode optional



URIT-8210

Automatic Chemistry Analyzer

- Random Access, 330T/H
- 90 reaction positions, 59 reagent positions, 71 sample positions, including detergent position, standard positions, QC positions, STAT positions
- High-accuracy grating optical system
- Durable ceramic syringes, 3 unattached probes, 2 high-performance mixers
- 24 hours non-stop water cycle cooling system to keep the reagent at 2-8°C
- Metal solid heating for the 37°C incubation
- Support LIS interface



URIT-880 Chemistry Analyzer

- Friendly operation with 7 inch colourful touch screen
- Sophisticated Software, large memory capacity
- High-quality filters & close optical system
- Convenient maintenance design
- 8 incubation positions



Point-of-care Testing Hemoglobin Meter

URIT-12 Hemoglobin Meter

- World's leading technology
- Fast hemoglobin and HCT test
- Accurate and reliable result
- Easy 3-step operation
- Suitable for various workplaces
- Hemoglobin control available



URIT-810 Chemistry Analyzer

- High-quality filter & close optical system
- Easy keypad operation
- Long-life halogen lamp
- Easy maintenance



URIT-15 Hemoglobin Meter

- Fast hemoglobin and HCT test within 10 seconds
- Accurate & reliable result
- Easy 3- step operation
- Hemoglobin control available



URIT-10

Uric Acid and Blood Glucose Meter

- High accuracy, only use 1.5ul whole blood
- Two units optional: $\mu\text{mol/L}$ and mg/dL for uric acid, mmol/L and mg/dL for blood glucose
- Automatic judgement for strips and code card



URIT-28 Glucose Meter

- Strip back slide button
- Easy and fast ,high-accuracy
- Glucose dehydrogenase method



URIT-25

Blood Glucose Meter

- Fast Blood Glucose test within 5 seconds
- High accuracy, only use 1ul whole blood
- Easy to use, automatic judgement for strips and code card



Urine Analyzer

URIT-31

Urine Analyzer

- Up to 120T/H
- Portable with touch screen
- Bluetooth printer and data transportation
- Up to 14 parameters with control solution
- Touch screen



Electrolyte Products

URIT-910C Plus Electrolyte Analyzer

- Test 6 parameters with 45 seconds: K, Na, Cl, Ca, pH, only use 100ul sample
- Individual reagent optional, cost-effective
- Friendly operation with touch screen
- Auto-sampling tray (25 positions) available
- Barcode scanner and bluetooth
- Intelligent bubble detection



MicroPlate Reader & Washer Products

URIT-660 Microplate Reader

- 8 channels, reading 96-well plate in 5 seconds
- Touch screen, large LCD display
- High accurate optical measuring system
- PC control software optional



Semi-coagulation Products

URIT-600/610 Coagulation Analyzer

- 2 channels(URIT-600)/4 channels (URIT-610) optional
- 16 sample positions, 4 reagent positions
- High-quality sample pipettor
- Advanced dual-magnetic circuit bead method
- Low reagent consumption



URIT-670 Microplate Washer



- Minimal residual volume (<2ul/well)
- 8-way and 12-way wash heads interchangeable
- Positive displacement pump with low noise
- Overfill washing and overflow protection
- Bottom washing function



Medical Software

LIS(Laboratory Information System)

- Easy installation & user-friendly interface
- Better patient care & more effective use of operation
- Convenient test result query & customized test report
- Improved quality of laboratory information management



Reagent

URIT Hematology Reagent & Control

- Available for: URIT hematology analyzer
- Available for: SYSMEX, ABBOTT, NIHON KOHDEN, ABX, MINDRAY, etc.
- **Hematology Controls: Bio-Rad control, URIT control**



URIT Urine Reagent Strip

- Available for Analyzers
- Available for Manual



URIT Chemistry Reagent

48 test items:

- TB/TB-VA/TB-DCA, DB/DB-VA/DB-DCA, TP, ALB, ALT, AST, GGT, ALP, CHE, TBA, PA, ADA
- Cr/CR-OX, Urea, UA, CSF/UTP, β 2-MG, MALB, Cys-c, CO₂
- CHOL, TG, HDL-C, LDL-C, Apo-B, HCRP, Lp (a)
- GLU-PAP/GLU-HK
- CK-NAC, α -HBDH, LDH, CK-MB, LDH1
- Fe, Ca, Mg, P
- ASO, RF, CRP
- IgG, IgA, IgM, C3, C4
- α -AMY, LPS, PAMY
- TF



URIT UQ Series Control


- Control for 10 parameters, includes level I and level II per kit
- Control for 11 parameters, includes level I, level II and level III per kit
- Control for 14 parameters, includes level I, level II and level III per kit



URIT Electrolyte Reagent

- Standard solution A, Standard solution B
- Internal filling solution, Reference solution
- Cleaning solution, Active solution
- QC solution, Protein-removing solution





Veterinary Product



URIT-5250Vet 5-Diff Hematology Analyzer

- 28 parameters including RET%, RET#, IRF
- 60T/H
- Laser light multi-dimensional cell classification
- Resistant RBC mode
- Support both whole blood and capillary blood samples
- Only three reagents for 25 parameters

URIT-3000Vet Plus Hematology Analyzer

- 16 pre-determined mode settings
- 60T/H, 21 parameters and 3 histograms
- 10.4 inch color LCD display & Linux operation system
- Support USB data backup & system upgrade
- Support LIS and HIS with HL7 protocol
- Touch screen optional



URIT-2900Vet Plus Hematology Analyzer

- 16 pre-determined mode settings
- 30T/H, 21 parameters and 3 histograms
- 10.4 inch color LCD display & Linux operation system
- Support USB data backup & system upgrade
- Support LIS and HIS with HL7 protocol
- Touch screen optional



URIT-50 Vet Urine Analyzer

- Up to 120T/H and 14 parameters available
- FDA approved
- Portable
- Barcode reader optional
- URIT UQ-10,UQ-11 and UQ-14 control available



URIT-880 Vet Chemistry Analyzer

- Friendly operation with 7 inch colourful touch screen
- Sophisticated Software, large memory capacity
- High-quality filters & close optical system
- Convenient maintenance design



URIT-31 Vet Urine Analyzer and Urine Strip

- Up to 120T/H
- Portable size with touch screen
- Bluetooth printer and data transportation
- Up to 14 parameters with control solution



URIT 11V Vet Urine Reagent Strip

- Easy visual reading
- Much higher reliability with Ascorbic Acid parameter



URIT-9000 Vet Integrated Diagnostic System

- General clinic information technology solution
- Wide angle detection
- Detachable disinfection
- Combination design
- Support bluetooth

Medical Electronic Product



URIT-9000 Integrated Diagnostic System (wall-mounted)

Blood Pressure, SpO2, Temp, Pulse, Otoscope, Ophthalmoscope,
Hammer, Nasal examination, Visual test, Color-blind test

- General clinic information technology solution
- Wide angle detection
- Detachable disinfection
- Combination design
- Support bluetooth



URIT-7200 Integrated Diagnostic System (wall-mounted)

Blood Pressure, SpO2, Temp, Pulse, Otoscope, Ophthalmoscope,
Hammer, Nasal examination, Visual test, Color-blind test

- Multifunction with 8 parameters
- Space-saving with integrated wall-mounted design
- Convenient for users with intelligent system
- More durable with innovative electronic display system