

BEST QUALITY

Since 2012



CE, ISO, GMP CERTIFIED

A-VAC TUBES











Our other products :

- ▶ Urine container - 30, 40, 50 & 60ml (Sterile / Non Sterile)
- ▶ Pipette tips - 200ml, 1000ml (Yellow & Blue)
- ▶ Cuvettes
- ▶ ESR Pippetes



A-VAC TUBES NON VACUUM

A-VAC PLUS VACUUM TUBES

Product	Photo	Blood Draw (ml)	Tube Size (mm)		Additive (s)	Specimen for Test	Laboratory Use
EDTA Tube		1, 2, 3, 4	12*75	PP	EDTA TAK3/K2	Whole Blood	For Hematology Determination, Blood Type Verification, Blood Cross-Match Test
Clot Activator		1, 2, 3, 4	12*75	PP	Clot Activator	Serum	Serum Biochemical Test
Sodium Citrate Tube 3.2		1.8	12*75	PP	3.2% Sodium Citrate	Plasma/ Whole Blood	For Blood Coagulation Determinations
ESR Tube 3.8		1.6	12*75	PP	3.8% Sodium Citrate	Whole Blood	For Corpuscle Sedimentation Test
Fluoride & Oxalate Tube		1, 2, 3, 4	12*75	PP	Potassium Oxalate & Sodium Fluoride	Plasma	For Blood Glucose Determinations
Plain Tube		1, 2, 3, 4	12*75	PP	No Additive	Serum Container	For Preservation of Blood Specimen
Lithium Heparin Tube		1, 2, 3, 4	12*75	PP	Gel & Lithium Heparin	Plasma	For Plasma Bio-chemical Determination, Hemorheology Test
Clot Activator		1, 2, 3, 4 & 5ml	12*75 12*100	PP	Gel & Clot Activator	Serum	For Serum Separation Biochemical Test

Non Vacuum Tubes 1 carton will be 12x100 No's

Vacuum Tubes 1 carton will be 18x100 No's



Our Products are being used by many customers from every part of the India

Manufactured by :

ACEE TUBES PVT. LTD.,

Panicker's Building, #14/21, School Street,
Sathya Nagar, Padi, Chennai-600 050.

Ph. : 044 2624 4545 Mob. : 90420 91711

e-mail : aceetubes@gmail.com / sales@aceetubes.com

www.aceetubes.com



Branches :

Tamilnadu & Andhra : 9042091712 | Kerala :9176941712 / 9176631712

Karnataka : 9176101712 | Pune : 9822218554 | Chandigarh : 904180820