




ELISA

A complete solution for diagnostic laboratory testing and blood screening



Your IVD
specialists!



A complete Panel of CE certified tests for the diagnosis of serological tropical diseases

DENGUE VIRUS

Dengue Virus IgG

Qualitative/semi-quantitative determination of IgG antibodies

Dengue Virus IgM

Qualitative determination of IgM antibodies

Dengue Virus NS1 Ag

Determination of Dengue Virus NS1 Antigen

TRYPANOSOMA CRUZI (Chagas)

T. Cruzi Ab

Third generation Enzyme Immunoassay for the determination of total antibodies to Trypanosoma Cruzi

WEST NILE VIRUS

West Nile Virus Antibodies

Qualitative determination of total antibodies

West Nile Virus IgG

Detection of IgG antibodies

West Nile Virus IgM

Detection of IgM antibodies

MALARIA

Malaria Ab

Third Generation Enzyme Immunoassay for the determination of total antibodies to Plasmodium species (falciparum, ovale and vivax)

CHIKUNGUNYA VIRUS

Chikungunya IgG (CHIKV IgG)

Detection of IgG antibodies

Chikungunya IgM (CHIKV IgM)

Qualitative determination of IgM antibodies

ZIKA VIRUS

ZIKA IgG (ZIKV IgG)

Qualitative determination of IgG antibodies

ZIKA IgM (ZIKV IgM)

Qualitative determination of IgM antibodies

ZIKA IgG Avidity Test (ZIKV IgG Avidity Test)

Determination of the avidity of IgG antibodies

YELLOW FEVER VIRUS

Yellow Fever Virus IgG

Detection of IgG antibodies for the follow-up of the vaccination

BLOOD SCREENING

The sole manufacturer of a **COMPLETE ELISA** kit panel
for Blood Bank screening requirements

HEPATITIS B VIRUS

HBs Ag one version ULTRA

Fourth generation Enzyme Immunoassay for the determination of Hepatitis B surface Antigen or HBs Ag (REC included)

HBc Ab

Determination of total antibodies to Hepatitis B core Antigen

HEPATITIS C VIRUS

HCV Ab

Version 4.0 Immunoassay for the determination of total antibodies

HTLV

HTLV I&II Ab Version ULTRA

Third generation Enzyme Immunoassay for the determination of antibodies to Human T-Cell Lymphotropic Virus type I&II

HIV

HIV 1,2 & 0 Ab&Ag (HIV Ab&Ag)

Fourth generation Enzyme Immunoassay for the determination of antibodies to HIV type 1&2&0 and P24 HIV-1 Antigen

MALARIA

Malaria Ab

Third Generation Enzyme Immunoassay for the determination of total antibodies of Plasmodium species (falciparum, ovale and vivax)

TRYPANOSOMA CRUZI (Chagas)

T. Cruzi Ab

Third Generation Immunoassay for the determination of total antibodies

TREPONEMA PALLIDUM (Syphilis)

Syphilis Ab

"Double sandwich" (ULTRA) immunoassay for the determination of total antibodies.

Available in one short step protocol

(Syphilis Ab One Version ULTRA)

or in two steps protocol optimized for automation

(Syphilis Ab Version ULTRA)

Also available...

DENGUE VIRUS

DENGUE Virus NS1 Ag

Determination of Dengue Virus NS1 Antigen

HEPATITIS E VIRUS

HEV Ab Version ULTRA

"Double sandwich" (ULTRA) immunoassay for the determination of total antibodies





A wide panel of CE certified Serology kits

HELICOBACTER PYLORI

HP IgG/IgM/IgA

Quantitative/qualitative determination of IgG/IgM/IgA antibodies

CagA IgA/IgG

Quantitative/qualitative determination of IgA/IgG antibodies to Helicobacter pylori Cytotoxin Associated Gene Antigen

HP Ag

Qualitative/quantitative determination of Helicobacter Pylori Antigen in human stools

CHLAMYDIA

Chlamydia Pneumoniae IgA/IgG/IgM

Determination of IgA/IgG/IgM antibodies to Chlamydia Pneumoniae

Chlamydia Trachomatis IgA/IgG/IgM

Determination of IgA/IgG/IgM antibodies to Chlamydia Trachomatis

PARVO VIRUS B19

Parvovirus B19 IgG/IgM

Quantitative/qualitative determination of IgG antibodies

MENINGITIS

Meningitis IgG

Determination of IgG antibodies

COXSACKIE B VIRUS

CoxB IgG

Determination of IgG antibodies

CoxB IgM

Determination of IgM antibodies

MEASLES VIRUS

Measles Virus IgG

Semi-quantitative determination of IgG antibodies

Measles Virus IgM

Qualitative determination of IgM antibodies

TREPONEMA PALLIDUM (Syphilis)

SYPH IgM

Determination of IgM antibodies

Syphilis Ab

"Double sandwich" (ULTRA) immunoassay for the determination of total antibodies.

Available in one short step protocol

(Syphilis Ab One Version ULTRA)

or in two steps protocol optimized for automation

(Syphilis Ab Version ULTRA)

HEPATITIS



A UNIQUE panel of CE certified kits for the diagnosis of Viral Hepatitis in human serum and plasma

HEPATITIS A

HAV Ab

Determination of total antibodies

HAV IgM

“Capture” Enzyme Immunoassay for the determination of IgM class antibodies

HEPATITIS B

HBc Ab

Determination of total antibodies to Hepatitis B core Antigen

HBc IgM

“Capture” Enzyme Immunoassay for the quantitative/qualitative determination of IgM class antibody to Hepatitis B Virus core Antigen

HBs Ab

Qualitative/quantitative determination of total antibodies to Hepatitis B surface Antigen

HBs Ag one version ULTRA

Fourth generation Enzyme Immunoassay for the determination of Hepatitis B surface Antigen

HBs Ag Confirmation

Set of Reagents for the confirmation of HBs Ag positivity

HBe Ag&Ab

Determination of Hepatitis B Virus “e” Antigen and Antibody

HEPATITIS C

HCV Ab

Version 4.0 Enzyme Immunoassay for the determination of total antibody

HCV IgM

Quantitative/qualitative determination of IgM antibodies

HCV Ab Confirmation

Microstrip based Enzyme Immunoassay for the confirmation of HCV Ab positivity

HEPATITIS D

HDV Ab

Competitive Immunoassay for the qualitative determination of total antibodies

HDV Ag

Determination of Hepatitis Delta Virus

HDV IgM

“Capture” Enzyme Immunoassay for the determination of IgM antibodies

HEPATITIS E

HEV Ab Version ULTRA

“Double sandwich” (ULTRA) immunoassay for the determination of total antibodies

HEV IgG

Determination of IgG antibodies

HEV IgM

Determination of IgM antibodies





A complete Panel of CE certified Prenatal Screening kits

CYTOMEGALOVIRUS

CMV IgG

Quantitative/qualitative determination of IgG antibodies

CMV IgM

“Capture” Enzyme Immunoassay
for the determination of IgM antibodies

CMV Antibody Avidity Test*

Determination of “Avidity” of IgG class antibodies

TOXOPLASMA GONDII

Toxo IgG

Quantitative/qualitative determination of IgG antibodies

Toxo IgM

“Capture” Enzyme Immunoassay
for the determination of IgM antibodies

Toxo Antibody Avidity Test*

Determination of “Avidity” of IgG class antibodies

RUBELLA VIRUS

RUB IgG

Quantitative/qualitative determination of IgG antibodies

RUB IgM

“Capture” Enzyme Immunoassay
for the determination of IgM antibodies

RUB Antibody Avidity Test*

Detection of the avidity of IgG antibodies

HERPEX SIMPLEX VIRUS

HSV1 IgG

Quantitative/qualitative determination of
IgG antibodies to Herpes Simplex Virus type 1

HSV1 IgM

“Capture” Enzyme Immunoassay for the determination
of IgM antibodies to Herpes Simplex Virus type 1

HSV2 IgG

Quantitative/qualitative determination of
IgG antibodies to Herpes Simplex Virus type 2

HSV2 IgM

“Capture” Enzyme Immunoassay for the determination
of IgM antibodies to Herpes Simplex Virus type 2

HSV1&2 IgG

Quantitative/qualitative determination of
IgG antibodies to Herpes Simplex Virus type 1&2

HSV1&2 IgM

“Capture” Enzyme Immunoassay for the determination
of IgM antibodies to Herpes Simplex Virus types 1&2

HUMAN PAPILLOMA VIRUS

Papilloma Virus IgG

Quantitative/qualitative determination
of IgG antibodies (R.U.O.)

PARVO VIRUS B19

Parvovirus B19 IgG

Quantitative/qualitative determination of IgG antibodies

Parvovirus B19 IgM

Qualitative determination of IgM antibodies

*NOT CE



A wide panel of CE certified Autoimmune disease kits

IgG anti SSB

Quantitative determination of IgG autoantibodies against SSB autoantigens

IgG anti SSA 60 KD

Quantitative determination of IgG autoantibodies against SSA-60KD autoantigens

IgG anti SSA 52 KD

Quantitative determination of IgG autoantibodies against SSA-52KD autoantigens

IgG anti U1-snRNP 68

Quantitative determination of IgG autoantibodies to RNP-68KDa

IgG anti dsDNA

Quantitative determination of IgG autoantibodies against double stranded DNA (dsDNA)

IgG anti Sm

Quantitative determination of IgG autoantibodies against Sm autoantigens

ANA Screening IgG

Qualitative determination of IgG antibodies to anti Nuclear Antigens

IgG anti Jo-1

Quantitative determination of IgG autoantibodies against Jo-1 autoantigens

ANA 8 parameters profile

Qualitative determination of IgG antibodies to Ro/SSA60KDa, Ro/SSA52KDa, La/SSB, RNP-68, Sm, Scl-70, CENP-B

IgG anti Centromere B

Quantitative determination of IgG autoantibodies against Centromere B autoantigens

ENA Screening IgG

Qualitative determination of IgG antibodies to Extractable anti Nuclear Antigens

IgG anti Scl-70

Quantitative determination of IgG autoantibodies against Scl-70 autoantigens

ENA 6 parameters profile

Qualitative determination of IgG antibodies against Ro/SSA, La/SSB, RNP-68, Sm, Scl-70, Jo-1

DIA.BLOod

ELISA Micro-plate Workstation

FULLY AUTOMATED HIGH THROUGHPUT



MAIN FEATURES

- ✓ DISPOSABLE TIPS (STANDARD)
- ✓ 9 INDEPENDENT INCUBATORS
- ✓ 3 INDEPENDENT WASHERS
- ✓ 8 CHANNEL INDEPENDENT MANIFOLD
- ✓ PRESSURE SENSOR & CLOT DETECTOR
- ✓ BAR-CODE IDENTIFICATION
- ✓ CHEMILUMINESCENCE READER (OPTIONAL)
- ✓ DISPENSING VERIFICATION SYSTEM (OPTIONAL)

ANALYZERS

- A) DIA.BLOod:
Optical Density Reader
- B) DIA.BLOod DUO:
Double Reading System Optical
Density & Relative Light Unit (CLIA)

CE IVD 98/79/EC
CB Certified

A CE certified combination of kits and instruments
For more information visit: www.diapro.it



Dia.Pro Diagnostic Bioprobes S.r.l.
Via Carducci, 27 - 20099 Sesto San Giovanni (MI) - Italy
Tel. +39 02 27.00.71.61/64.50 Fax +39 02 44386771
<http://www.diapro.it> - email: info@diapro.it

EN ISO 13485 Certified

BRO-ELISA Rev. 1018