







COVID-19

IgG/IgA/IgM EIA

-  Qualitative detection of COVID (SARS-Cov-2) IgG/IgA/IgM EIA Antibody
-  Principle based on an indirect solid-phase enzyme immunoassay
-  Sample Type – Serum / Plasma
-  Available in convenient pack size of 96 Tests
-  Longer Shelf Life: 12 months at 2° C - 8° C



Flexibility in Operation

- The Kit can be used for diagnostic as well as serological surveillance purposes

Certified Quality

- CE IVD marked reliable quality

Additional Information

- Compatible with major semi & fully automated ELISA analysers
- The Kit is based on an indirect solid-phase enzyme immunoassay with horseradish peroxidase as a marker enzyme.
- Colour coded reagents for ease of identification
- Negative & Positive Control included in kits

Kit Content

- COVID-19 Antigen Coated Microplate
- Sample Diluent
- Negative Control
- Cut-Off Control
- Positive Control
- Conjugate
- TMB / Substrate Solution
- Stopping Solution
- Washing Solution

Availability

- Immediate in convenient kit size of 96 tests

Product Code:

Covid IgG EIA	6111700
Covid IgA EIA	6111710
Covid IgM EIA	6111720

**For more information please order your IFU from
Labsystems Diagnostics Oy**

Tiilitie 3, FI-01720 VANTAA, Finland Tel: +358 (0) 20 155 7530
sales@labsystemsdx.com | www.labsystemsdx.com